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DONG Hengtao, AI Yun, LI Yueqi, et al. Determination of sulfonamides in water from aquaculture environment by UHPLC-MS/MS[J]. Environmental Chemistry, 2023, 42 (2): 671-674.

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### Determination of sulfonamides in water from aquaculture environment by UHPLC-MS/MS

DONG Hengtao<sup>1</sup> AI Yun<sup>2</sup> LI Yueqi<sup>1</sup> DENG Xiaoli<sup>1</sup> PENG Yifang<sup>1</sup>  
HUANG Taohong<sup>1</sup>

(1. Shimadzu (China) Co., Ltd., Xi'an, 710065, China; 2. Xi'an Institute for Food and Drug Control, Xi'an, 710054, China)

**Abstract** To establish a quick and sensitive analytical method for the determination of 21 sulfonamides in water samples from aquaculture environment by ultra high performance liquid chromatography-electrospray tandem mass spectrometry (UHPLC-MS/MS). The water samples were extracted and cleaned by Solid Phase Extraction through Oasis HLB columns, and then Velox SP-C<sub>18</sub> column was used for separation. The target compounds were determined by multiple reaction monitoring (MRM) mode and quantified by internal standard method. It was shown that the analysis of one sample needed only 12.5 minutes. The calibration curve of 21 sulfonamides residues were constructed with correlation coefficients (*r*) more than 0.995, respectively. Good reproducibility on both retention time and peak area were observed. The quantitation limit of 21 sulfonamides were all 0.63 ng·mL<sup>-1</sup> in this method. The method could be used for fast and accurate quantitative determination of 21 sulfonamides residues in water.

**Keywords** tandem mass spectrometry, water, sulfonamides, solid phase extraction.

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### 参考文献 (References)

- [1] LI L, LIU D, ZHANG Q, et al. Occurrence and ecological risk assessment of selected antibiotics in the freshwater lakes along the middle and lower reaches of Yangtze River Basin [J]. *Journal of Environmental Management*, 2019, 249: 1-8.
- [2] LI W H, GAO L H, SHI Y L, et al. Occurrence, distribution and risks of antibiotics in urban surface water in Beijing, China [J]. *Environmental Science*, 2015, 17(9): 1611-1619.
- [3] 唐娜, 张圣虎, 陈政宏, 等. 长江南京段表层水体中 12 种磺胺类抗生素的污染水平及风险评价 [J]. *环境化学*, 2018, 37(3): 505-512T.
- [4] 李丹, 沈昕, 张玉洁, 等. 高效液相色谱法测定鸡肉中 5 种磺胺类药物及乙酰胺苯甲酯增效剂残留量 [J]. *食品安全质量检测学报*, 2020, 11(1): 202-208.
- [5] 李柳毅, 范辉, 范磊, 等. 固相萃取-高效液相色谱法测定地表水中 4 种磺胺类抗生素 [J]. *化学分析计量*, 2017, 26(6): 38-40.
- [6] XU Y, DING J, CHEN H Y, et al. Fast determination of sulfonamides from egg samples using magnetic multiwalled carbon nanotubes as adsorbents followed by liquid chromatography-tandem mass spectrometry [J]. *Food Chemistry*, 2013, 140: 83-90.
- [7] CUI H Y, CHANG H, ZHENG H J, et al. Determination and occurrence of sulfonamide transformation products in surface waters [J]. *Science of the Total Environment*, 2021, 799: 146562.
- [8] CHEN Y P, XIA S D, HAN X Q, et al. Simultaneous determination of malachite green, chloramphenicols, sulfonamides, and fluoroquinolones residues in fish by liquid chromatography-mass spectrometry [J]. *Journal of Analytical Methods in Chemistry*, 2020: 3725618.
- [9] 陈曼娜. 磺胺类药物的多残留酶联免疫法检测 [J]. *质量安全与检验检测*, 2020, 30(4): 46-47.
- [10] 中华人民共和国农业农村部. SC/T 9436-2020 水产养殖环境(水体、底泥)中磺胺类药物的测定[S]. 北京: 中国农业出版社, 2020.